



Surface Characterization of **Medical Devices**
Bridging metrology and healthcare

SENSOFAR

INTRODUCTION	3
--------------	---

CHALLENGES OF THIS INDUSTRY	4
-----------------------------	---

HOW CAN WE HELP	5
-----------------	---

Polivalent equipment	5
----------------------	---

Equipment you can trust	6
-------------------------	---

Not just suppliers, but partners	8
----------------------------------	---

METROLOGY IN PROSTHETICS	10
--------------------------	----

Hip implants	11
--------------	----

Knee implants	17
---------------	----

Dental implants	20
-----------------	----

Wedge implants for feet	23
-------------------------	----

CHARACTERIZATION OF OFTALMIC LENSES	27
-------------------------------------	----

Specular lenses for children	28
------------------------------	----

SURGICAL EQUIPMENT	30
--------------------	----

Intraoperation devices	31
------------------------	----

CARDIOVASCULAR DEVICES INSPECTION	33
-----------------------------------	----

Coronary stents	34
-----------------	----

Heart valves	35
--------------	----

Drug-eluting devices	36
----------------------	----

KEYS TO SENSOFAR'S SUCCESS	37
----------------------------	----





HOW CAN WE HELP

Polivalent equipment


Our systems are designed with inherent versatility, adapting easily to changing measurement requirements through advanced technology and numerous automation tools. This flexibility maximizes our equipment's long-term value.

Multiple measuring techniques in one system

Most of our systems are equipped with multiple measuring technologies, providing a comprehensive characterization and enabling users to select the best technology based on their application needs.

Easy transformation from R&D to QC

Our systems offer tremendous flexibility by easily transitioning from research and development (R&D) to quality control (QC). By utilizing specific features, our systems can meet the needs of both R&D and QC, providing users with a versatile solution for various applications.



F-operator: Leveling plane Algorithm: Gaussian
 S-filter λ_s : 0.0008333 mm L-filter λ_c : 0.25 mm


ISO 25178 / Functional stratified surface

Sk	1.9053 μm
Spk	0.8020 μm
Svk	0.1993 μm

Roughness for durability


During an implant's service life, wear leads to the formation of debris that can harm patient health by causing osteolysis and inflammatory reactions. The wear behavior of hip implants, influenced by the human gait cycle, is closely linked to their surface condition. Reducing wear and debris is crucial, particularly for younger, active patients with longer life expectancies. Research¹ has shown that measuring the roughness of the top area of acetabular cups evaluates their wear characteristics and, thus, their durability.

Confocal technology, ideal for imaging local slopes and dense features, was employed at a 20X magnification. Studies identified three key roughness parameters—core roughness depth (Sk), reduced valley height (Svk), and reduced peak height (Spk)—that correlate machining parameters with wear behavior, indicating that higher values lead to greater wear volumes.







When measuring roughness on a tilted sample, the results can be useful if used qualitatively—for example, when comparing different implants measured in the same region. However, these measurements are not quantitatively accurate for the following reasons.

- 1 The effective measured area changes with the tilt; for example, at a 30-degree tilt, the imaged area is 13% smaller than at 0 degrees.



- 2 The inclination of a sample significantly affects the measured roughness, as the sampling changes with different inclinations.





Coronary stents

Visual inspection capabilities allow for detecting defects and irregularities, which could compromise the device's performance or safety.

The [Q vix](#) system autonomously gathers information for the dimensional measurements and visual inspection of stents, enhancing efficiency and precision in the production process. This advanced system captures high-resolution images to meticulously inspect both the outer and inner surfaces of each stent, promptly identifying any irregularities. By operating autonomously, the [Q vix](#) system significantly accelerates sample throughput on the production line, ensuring consistent quality control and reducing inspection time.

Unlock the complete
White paper
by signing up to

my SENSOFAR

